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1. Lillie, R.D. and Ashburn, L.L. 1943. Super-saturated solutions of fat stains in dilute isopropanol for demonstration of acute fatty degenerations not shown by Herxheimer technique. Arch. Pathol. 36: 432.
2. Moscatello, D.K. Cell culture media, kits, and methods of use. U.S. Patent Application 11542863, filed 10/04/2006; U.S. Provisional Patent Application No. 60/723,804, filed 10/06/2005.
3. Pittenger M.F., Mackay A.M., Beck S.C., Jaiswal R.K., Douglas R., Mosca J.D., Moorman M.A., Simonetti D.W., Craig S., and Marshak D.R. 1999. Multilineage potential of adult human mesenchymal stem cells. Science 284(5411):143-147.
4. Sekiya, I., Larson, B.L., Smith, J.R., Pochampally, R., Cui, J-G., and Prockop, D.J. 2002. Expansion of Human Adult Stem Cells from Bone Marrow Stroma: Conditions that maximize the yields of early progenitors and evaluate their quality. Stem Cells 20:530-541.
5. Sheehan, D.C. and Hrapchak, B.B. 1980. Lipids. In Theory and Practice of Histotechnology, Second Edition, Sheehan, D.C. and Hrapchak, B.B., Eds., Battelle Press, Columbus, OH. Pp202213.

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1. Alizadeh, Effat, et al. "Upregulation of MiR 122 via trichostatin a treatments in hepatocyte like cells derived from mesenchymal stem cells." Chemical biology & drug design 87.2 (2016): 296-305.
2. Abolhassani, S., K. Yavari, and J. Mohammadnejad. "Investigation of Isolating of Human Mesenchymal Stem Cells from Umbilical Cord Blood and Labeling them with 99m Technetium." Int J Stem Cell Res Transplant 4.6 (2016): 190-194.
3. Jalilzadeh-Tabrizi, Sepideh, et al. "A Biomimetic Emu Oil-Blended Electrospun Nanofibrous Mat for Maintaining Stemness of Adipose Tissue-Derived Stem Cells." Biopreservation and biobanking 16.2 (2018): 66-76.
4. Fayazi, Mehri, Mojdeh Salehnia, and Saeideh Ziaei. "Differentiation of human CD146-positive endometrial stem cells to adipogenic-, osteogenic-, neural progenitor-, and glial-like cells." In Vitro Cellular & Developmental Biology-Animal 51.4 (2015): 408-414.
5. Davoodian, Nahid, et al. "MicroRNA 122 overexpression promotes hepatic differentiation of human adipose tissue derived stem cells." Journal of cellular biochemistry 115.9 (2014): 1582-1593.
6. Alizadeh, Effat, et al. "The effect of dimethyl sulfoxide on hepatic differentiation of mesenchymal stem cells." Artificial cells, nanomedicine, and biotechnology 44.1 (2016): 157-164.
7. Dayer, Dian, et al. "Sonic hedgehog pathway suppression and reactivation accelerates differentiation of rat adipose-derived mesenchymal stromal cells toward insulin-producing cells." Cytotherapy 19.8 (2017): 937-946.
8. Taghikani, Mohammad, and Fatemeh Eskandari. "Differentiation of Human MSCs Into Nsulin Producing Cells by Using Lentiviral Vector Carrying PDX-1."